

QtDaq python analysis tools

Chloé Sole

Updated: April 5, 2022

This document covers the base workings of the `read_dat` class and `event` class developed for analysis of `dat` files produced by QtDaq software [ref].

To use these classes please ensure the associated python files are in the same folder as your script and that the dependencies are installed.

Contents

1	Class <code>read_dat</code>	2
1.1	Parameters	2
1.2	Attributes	2
1.3	Public Methods	3
1.3.1	<code>read_event</code>	3
1.3.2	<code>lst_out</code>	3
1.3.3	<code>get_fails</code>	4
1.3.4	<code>add_selections</code>	4
1.3.5	<code>select_events</code>	5
1.4	Private Methods	5
1.4.1	<code>onclick</code>	5
1.4.2	<code>press</code>	6
1.5	Fails	6
2	Class <code>event</code>	6
2.1	Parameters	7
2.2	Attributes	7
2.3	Public Methods	7
2.3.1	<code>get_event_id</code>	7
2.3.2	<code>get_ch</code>	8
2.3.3	<code>get_t0</code>	8
2.3.4	<code>get_trace</code>	8
2.3.5	<code>get_CFD</code>	8
2.3.6	<code>get_baseline</code>	8
2.3.7	<code>get_long_integral</code>	8
2.3.8	<code>get_short_integral</code>	9
2.3.9	<code>get_pulse_shape</code>	9
2.3.10	<code>get_pulse_height</code>	9
2.3.11	<code>get_fails</code>	9
2.4	Private Methods	10
2.4.1	<code>cfid</code>	10
2.4.2	<code>sum_integral</code>	10
2.4.3	<code>check_polarity</code>	10
2.5	Fails	10
3	Examples	11

1 Class read_dat

Parameters: file_name, sample_rate, CFD, t_start, t_long, t_short, baseline_samples, output

Dependencies: scipy, event, matplotlib, csv

Public Attributes: headerSize, maxChannels, preambleSize, fileName, inputFile, header, eventCounter, eventTimeStamp, endFile, nsPerSample, CFD (array of length 3), tStart, tShort, tLong, chActive, fails, totFails, baselineSamples, selection, cuts

Public Methods: read_event, lst_out, get_fails, add_selections, select_events

Private Methods: onclick, press

This class provides an object representing a single file. The file is assumed to have an event structure as follows:

headerSize: 72 bytes
maxChannels: 64 ch
preambleSize: $4+20+4*\text{maxChannels}$ bytes

Each event is preceded by a preamble containing an event timestamp ($\text{preamble}[5]*8\times 10^3 \mu\text{s}$) and the size of each channel's acquisition ($\text{preamble}[6:]$).

The read_dat class needs to be imported into your script in order to use these tools.

```
from read_dat import *
```

To then create an object call read_dat()

```
file = read_dat("filename.dat",...)
```

1.1 Parameters

read_dat(*self*, file_name, sample_rate, CFD, t_start, t_long, t_short, baseline_samples, output)

file_name: *str* of the file name must include file suffix ".dat"

sample_rate: float, ns per sample for the acquired data.
Default: 2 ns (DT5730 sample rate)

CFD: 1D array of length 3: [Fraction, Length (samples), Offset (samples)]
Refer to event.CFD for more information on what these parameters are
default: [0.75, 6, 6]

t_start: start time (ns) from the CFD defined t_0
default: -80 ns

t_long: long integral gate (ns) from the CFD defined t_0
default: 400 ns

t_short: short integral gate (ns) from the CFD defined t_0
default: 10 ns

baseline_samples: length of baseline calculation in samples
default: 200 samples

1.2 Attributes

All attributes are currently public. Please edit attributes directly with **caution** as doing this will not update events already processed in the file.

headerSize: *int*, header size of 72 bytes
maxChannels: *int*, maximum number of channels allowed 64
preambleSize: *int*, $4+20+4*\text{maxChannels}$ bytes
fileName: *str*, input file name
inputFile: *file*, the open binary file
header: *binary*, header for the open file

<code>eventCounter:</code>	<i>int</i> , running event counter to keep track of where you are in the file
<code>eventTimeStamp:</code>	<i>float</i> , time stamp in μ s of the last event read
<code>endFile:</code>	<i>boolean</i> , if end of file has been reached <code>endFile = True</code>
<code>nsPerSample:</code>	<i>float</i> , the time in ns between samples
<code>CFD:</code>	<i>1D array of length 3</i> , parameters required to determine a Constant Fraction Discriminator filter of the events [Fraction, Length (samples), Offset (samples)]
<code>tStart:</code>	<i>int</i> , start time (ns) from the CFD defined t_0
<code>tShort:</code>	<i>int</i> , short gate end time (ns) from the CFD defined t_0
<code>tLong:</code>	<i>int</i> , long gate end time (ns) from the CFD defined t_0
<code>chActive:</code>	<i>1D array</i> , an array of the active channels indices
<code>fails:</code>	<i>$n \times 5D$ array</i> , fail tracker [start, long, short, integrals, zero] per ch, a description of the fail checks can be found in section 1.5.
<code>totFails:</code>	<i>1D array</i> , running count of the number of failed events per channel
<code>baselineSamples:</code>	<i>int</i> , number of samples used to calculate the baseline, taken from the first sample in the event
<code>selection:</code>	<i>2D list</i> , list of x and y co-ordinates used for enclosed area selections
<code>cuts:</code>	<i>1D array</i> , <code>len(cuts)</code> is equal to the number of selections available, with each element containing the number of co-ordinates in that selection

1.3 Public Methods

1.3.1 read_event

`read_dat.read_event(self)`

Method to read one event from the open file.

Parameters:	No additional parameters
Returns:	ev: <i>event</i> array Array of <i>event</i> objects associated with that event just read for the active channels. Returns <i>True</i> when the file end is reached
Raises:	

1.3.2 lst_out

`read_dat.lst_out(self, events=False, ch=True, output=True, traces=False, cuts=False, filename=)`

Method to read multiple events or whole file from the open file with list mode output options. Produces a **params** and **trace** csv file per channel if requested.

Parameters:	events: <i>int</i> or <i>boolean</i> , optional The number of events to read and output the desired list mode information. If <i>False</i> the full file is read. <i>Default</i> value is <i>False</i> .
	ch: <i>int</i> array or <i>boolean</i> , optional If <i>True</i> all channels are read out. If <i>int</i> array only the selected channels are read out (channel numbering from 0). <i>Default</i> value is <i>True</i> .
	output: <i>int</i> array or <i>boolean</i> , optional If <i>True</i> L [ch], S [ch], T_{trigger} [us], baseline [bits], pulse height [bits] are read out into a file per channel. If a binary <i>int</i> array (len=5) only the selected parameters are read out. The format of the array is [L,S,T,baseline,PH] where 1 is an indication to output that parameter and 0 is to ignore that parameter. <i>Default</i> value is <i>True</i> .
	traces: <i>boolean</i> , optional If <i>False</i> no traces will be outputted.

If **True** traces will be outputted list mode in a csv file per channel.

Default value is **False**.

cuts: int array or boolean, optional

If **False** no cuts are taken into account in the outputted data.

If an **int** array (len = number of cuts added) where the index of the array is the cut id and the value being 1, 0 or -1 includes, ignores or excludes that cut.

Default value is **False**.

filename: str, optional

Desired output file name.

Default value is empty which produces *fileName__params__ch.csv*.

Output: *fileName__params__ch.csv*: a csv file per channel with the list mode output of the analysed parameters for the desired number of events.

fileName__traces__ch.csv: a csv file per channel with the list mode output of the traces (in bits) for the desired number of events.

Raises: No additional error codes.

1.3.3 get_fails

`read_dat.get_fails(self, display=False)`

Returns fail information for the processed events.

Parameters: **display**: boolean, optional

If **False**, nothing is printed out in the **cmd**,

If **True**, a summary of the fails per channel is printed to **cmd**

Default value is **False**.

Returns: **fails**: $n \times 5D$ int array for n ch,

in the format of [start, long, short, integral, zero] for n channels

The value at the associated index indicates the number of fails out of the processed events that have failed that check. Refer to section 1.5 for a breakdown of the fail checks.

totFails: int array, length = number of active channels

Number of events failed per channel

eventCounter: int

Total number of events processed

Output: If **display** is **True**, writes the break down of fails per channel to **cmd**.

Raises:

1.3.4 add_selections

`read_dat.add_selections(self, L=[], S=[], mode="m", lims = [[0, 50000], [0, 1]], file=False)`

Method to add multiple cuts to the events. Can be run in manual **m** or predetermined **p** modes. Requires lists of **L** and **S** values in mode **m** for visual support to determine the selections. No lists of **L** and **S** values are required for mode **p** and no visual will be displayed.

Parameters: **L**: float array, optional

Array of *L* values for the processed events

S: float array, optional

Array of *S* values for the processed events

mode: char, optional

m: Manual mode, includes visual assistance when setting the enclosed area selections, refer to the key presses in table 1.

p: Predetermined mode, allows one to enter past manual selections without the visual assistance. The format of this file should match the output file as described below.

Default value is **m**.

lims: 2x2 float array, optional

Array of x and y limits for the visual aid SL plot.

Default value is `[[0,50000],[0,1]]`.

file: False, or *str*, optional

if **false** no output file is set for the selections and they are not saved unless.

if mode is **p** then a file name is required *str*, this is taken as the file name for the input file for the predetermined selections.

Default value is **False**.

Output: If **file** is a *str* and mode is **p** then an output file is created *fileName__cuts__SL.csv*. The output file has a header of one line, the cuts are then written to the rows in pairs. The first row is the x co-ordinates for the first cut and the second row in the y co-ordinates for the first cut.

Returns: Returns nothing, but sets the attributes **selection** and **cuts** for the file.

Raises: No additional errors.

1.3.5 select_events

```
read_dat.select_events(self,L,S, cut_id=[0],inc=[1],visual=False,lims = [[0,50000],[0,1]])
```

Method to pull the events which fall within the desired area. This area is defined by the inclusion or exclusion of the added selections. **If no selections have been added to the read_dat object all events will be returned.** There is a visual aid to ensure you are selecting for the correct events.

Parameters: **L:** float array, required

Array of *L* values for the processed events

S: float array, required

Array of *S* values for the processed events

cut_id: int array, optional

array of cut id's which are either included or excluded.

Default value is `[0]`, referring to cut 0.

inc: int array, optional

array of 1's or -1's to indicate which cut is either included or excluded respectively, `len(inc) = len(cut_id)`.

Default value is `[1]`, including cut 0 as default.

visual: boolean, optional

If **True** a SL plot with the included and excluded events is displayed with the cut boundaries.

Default value is **False**.

lims: 2x2 float array, optional

Array of x and y limits for the visual aid SL plot.

Default value is `[[0,50000],[0,1]]`.

Returns: Returns *S* and *L* values which fall within the defined area.

Raises: No additional errors. **add check for validity of parameters**

1.4 Private Methods

1.4.1 onclick

```
read_dat.__onclick(self)
```

Method used by `read_dat.add_selections` through `read_dat.__press`. Determines the clicked co-ordinates and stores them in attribute **selection**. Accessible within the `read_dat` class. Is called by a key press **a** or **A** while in `read_dat.add_selections`.

1.4.2 press

`read_dat.__press(self)`

While in manual mode of `read_dat.add_selections`, the commands for adding, editing or removing selections are given by the key presses in table 1. In an active selection, **a** or **A** has been pressed, use the left mouse click to select co-ordinates. The selected co-ordinates will be printed to screen and saved to attribute `selection` of the associated file.

Table 1: The relevant key presses while in manual mode for adding selections

Key	Action
a, A	Add a new selection
u, U	Undo previous point, only usable while in a selection
x, X	End current selection, can only end a selection if there are more than 2 co-ordinates in the selection
d, D	Delete previous completed selection
q, Q	Quit, ends visual guide and re-enters the main code segment
o, O	Outputs the selections added to a file (<i>fileName_cuts_SL.csv</i>)

1.5 Fails

The array returned by `read_dat.get_fails` is a $n \times 5$ D array for n channels. The five fail conditions presented per channel are indicated in table 2.

Table 2: Description of the fail checks returned by `read_dat.get_fails`.

Index	Fail Name	Fail Condition
0	start	The start time is set outside of the acquisition window
1	long	The long integral end gate is outside of the acquisition window
2	short	The short integral end gate is outside of the acquisition window
3	integral	The calculated short integral is negative or the calculated long integral has a smaller value than the calculated short integral
4	zero	The CFD calculation failed to return a reasonable t_0

2 Class event

Parameters: <code>event_id</code> , <code>ch</code> , <code>t0</code> , <code>trace</code> , <code>CFD</code> , <code>integrals</code> , <code>baseline</code>
Dependencies: <code>numpy</code>
Public Attributes: <code>eventID</code> , <code>ch</code> , <code>triggerTime</code> , <code>baseline</code> , <code>trace</code> , <code>CFD</code> , <code>istart</code> , <code>ishort</code> , <code>ilong</code> , <code>longIntegral</code> , <code>shortIntegral</code> , <code>fails</code>
Public Methods: <code>get_event_id</code> , <code>get_ch</code> , <code>get_t0</code> , <code>get_trace</code> , <code>get_CFD</code> , <code>get_baseline</code> , <code>get_long_integral</code> , <code>get_short_integral</code> , <code>get_pulse_shape</code> , <code>get_pulse_height</code> , <code>get_fails</code>
Private Methods: <code>cfd</code> , <code>sum_integral</code> , <code>check_polarity</code>

This class provides an object representing a single event.

The `event` class needs to be imported into your script in order to manipulate event objects directly and not only through the `read_dat` class. If you are using the `read_dat` class you do not need to import the `event` class additionally as it is inherited. To import the `event` class use the following command:

```
from event import *
```

To then create an event object call `event()`, refer to section 2.1 for a description of the parameters you can use to initiate an `event` object.

```
ev = event(...)
```

2.1 Parameters

`event(self, event_id, ch, t0, trace, CFD, integrals, baseline)`

<code>event_id:</code>	<i>int</i> , event counter across all channels
<code>ch:</code>	<i>int</i> , which channel the event occurred on
<code>t0:</code>	<i>float</i> , trigger time in μ s
<code>trace:</code>	<i>1D array</i> , full y values of the trace in bits, the x values are the trace can be determined through the properties of the digitiser used to acquire the data
<code>CFD:</code>	<i>1D array</i> , an array of length 3 containing the CFD analysis parameters required: [fraction,length,offset]. Refer to the private method <code>cfid</code> (sec. []) for a break down of these parameters.
<code>integrals:</code>	<i>1D array</i> of length 3, containing the indices for the start and ends of the integrals in the format of [start, short, long]
<code>baseline:</code>	<i>int</i> , the number of samples used for the baseline calculation

2.2 Attributes

All attributes are currently public.

<code>eventID:</code>	<i>int</i> , event ID
<code>ch:</code>	<i>int</i> , channel the event was acquired on
<code>triggerTime:</code>	<i>float</i> , the trigger time in μ s
<code>baseline:</code>	<i>float</i> , the baseline as calculated as the average of the number of samples from the start of the acquisition window as provided by the user
<code>trace:</code>	<i>1D array</i> , array representing the <i>y</i> -values in bits for the sampled trace
<code>CFD:</code>	[<i>1D array</i> , <i>float</i>], the array represents the event post CFD filter and the returned float (<code>izero</code>) is the determined weighted average for the index of the zero-crossing
<code>istart:</code>	<i>int array</i> , the start index for the integrals, determined as <code>istart = izero + integrals[0]</code>
<code>ishort:</code>	<i>int array</i> , the end index for the short integral, defined as <code>ishort = izero + integrals[1]</code>
<code>ilong:</code>	<i>int array</i> , the end index for the long integral, defined as <code>ilong = izero + integrals[2]</code>
<code>longIntegral:</code> within	<i>float array</i> , the value of the long integral, determined as the sum of samples the user defined gates
<code>shortIntegral:</code> within	<i>float array</i> , the value of the short integral, determined as the sum of samples the user defined gates
<code>fails:</code>	<i>1D array</i> , <i>len</i> = 5, fail tracker [start, long, short, integrals, zero], a description of the fail conditions can be found in section 2.5

2.3 Public Methods

2.3.1 `get_event_id`

`event.get_event_id(self)`

Method to return the associated event ID.

Parameters:	none
Returns:	eventID: int The event ID of the event used to call this method.

2.3.2 `get_ch`

`event.get_ch(self)`

Method to return the channel on which the associated event was acquired.

Parameters: none

Returns: **ch**: int

The channel on which the event, used to call this method, was acquired.

2.3.3 `get_t0`

`event.get_t0(self)`

Method to return the trigger time in μs of the associated event.

Parameters: none

Returns: **triggerTime**: float

The trigger time (μs) of the event used to call this method.

2.3.4 `get_trace`

`event.get_trace(self)`

Method to returns the y -values of the trace for the associated event.

Parameters: none

Returns: **trace**: 1D int array

The y -values of the trace of the event.

2.3.5 `get_CFD`

`event.get_CFD(self)`

Method to return the CFD of the associated event.

Parameters: none

Returns: **cfdArr**: 1D float array,

The array representing the CFD filtered trace of the event.

zero_cross: float,

The index value of the zero-crossing, determined as a weighted average of the indices on opposite sides of the zero-point.

2.3.6 `get_baseline`

`event.get_baseline(self)`

Method to return the baseline of the associated event.

Parameters: none

Returns: **baseline**: float

Average of the first x , user defined, many samples of the trace.

2.3.7 `get_long_integral`

`event.get_long_integral(self)`

Method to return the long integral or integrals of the associated event.

Parameters: none

Returns: **longIntegral**: int

Sum of the samples within the user defined start and end gates for the long integral.

2.3.8 `get_short_integral`

`event.get_short_integral(self)`

Method to return the short integral or integrals of the associated event.

Parameters: none

Returns: **shortIntegral**: int

Sum of the samples within the user defined start and end gates for the short integral.

2.3.9 `get_pulse_shape`

`event.get_pulse_shape(self)`

Method to return the pulse shape value or values of the associated event determined through a charge comparison (CC) calculation. If there are multiple short and long gates pulse shape cannot be calculated.

The pulse shape parameter S , as determined through a CC, is defined as:

$$S = \frac{Q_s}{Q_l} \quad (1)$$

where, $Q_{s/l}$ represent the short and long integrals respectively.

Parameters: none

Returns: **S**: float

The pulse shape parameter S of the event as determined through the desired method.

2.3.10 `get_pulse_height`

`event.get_pulse_height(self)`

Method to return the pulse height of the associated event.

Parameters: none

Returns: **eventID**: float

The maximum of the measured trace.

2.3.11 `get_fails`

`event.get_fails(self, display=False)`

Returns fail information for the associated event.

Parameters: **display**: boolean, optional

If **False**, nothing is printed out in the **cmd**,

If **True**, a summary of the fails for the selected event is printed to **cmd**

Default value is **False**.

Returns: **fails**: $n \times 5$ D int array for n ch,

in the format of [start, long, short, integral, zero] for n channels

The value at the associated index indicates the number of fails out of the processed events that have failed that check. Refer to section 1.5 for a breakdown of the fail checks.

Output: If **display** is **True**, writes the break down of fails per channel to **cmd**.

2.4 Private Methods

2.4.1 cfd

`event.__cfd(self, F, L, O)`

Method used when an event is created to assign values to the attribute **CFD**. This filter is defined in equation 2 where v_i is the CFD filtered event sample i , L is the filter length in samples, F is the filter fraction, O is the filter offset and V is the voltage of the raw unfiltered signal, N and O are in either ns or samples as long as the units used are consistent throughout the calculation.

$$v_i = \sum_{j=1}^L (F V_{i-j} - V_{i-j-O}) \quad (2)$$

Parameters: **F**: float, required
 The fraction F which is used to scale the trace in equation 2.
 L: int, required
 The length of the filter L
 O: int required
 The offset of the filter O
 Returns: **cfdArr**: 1D float array,
 The trace array post CFD filter.
 zero_cross: float,
 The weighted sum of the indices of the samples closest to the zero-crossing.

2.4.2 sum_integral

`event.__sum_integral(self, end)`

Method used when an event is created to assign values to the long and short integrals.

Parameters: **end**: int, required
 The end index for the sum integral.
 Returns: **sum_int**: int,
 The sum of the trace from the start index up to and not including the end index.

2.4.3 check_polarity

`event.__check_polarity(self)`

Method used when initiating an event to determine the polarity of the pulse and to switch it if necessary for consistency in the calculations.

2.5 Fails

The array returned by `event.get_fails` is a 1×5 D array for a single event. The five fail conditions presented for the calculations a single event can fail are indicated in table 3.

Table 3: Description of the fail checks returned by `event.get_fails`.

Index	Fail Name	Fail Condition
0	start	The start time is set outside of the acquisition window
1	long	The long integral end gate is outside of the acquisition window
2	short	The short integral end gate is outside of the acquisition window
3	integral	The calculated short integral is negative or the calculated long integral has a smaller value than the calculated short integral
4	zero	The CFD calculation failed to return a reasonable t_0

3 Examples